

Digital Radiography for the Mobile World.



ViZion™ + DR Mobile

The ViZion™+ DR mobile kit features the lightweight wireless cesium panel along with the reliable SR130 portable X-ray system. The digital package uses a mounted laptop tray and convenient panel slide to hold the laptop and the panel.

ViZion™ + Wireless DR Specifications

Weight: 7.9 lbs (CSI)
Carbon Fiber and Mg Alloy Construction
Pixel Matrix - 2304 x 2800
3.7 lp/mm (CSI)
Cassette-sized
Lithium Ion Battery (2 Included)

SR 130 X-ray System Specifications

3.0kw Output - Fully Regulated mA / kVp
KVP Range 40-100
MA Range 15 and 30
MAS Range 0.15-120

ViZion is a trademark of Konica Minolta Healthcare Americas, Inc.
© 2017 Konica Minolta Healthcare Americas, Inc.