



# Imaging Services

DR & CR DIGITAL EQUIPMENT

X-RAY SUPPLIES AND SERVICE

## SPECIAL DISCOUNTED PRICING

# RADIATION ATTENUATING GLOVES

Our radiation attenuating surgical gloves are designed to protect the physician or surgeon's hands while allowing for maximum flexibility, dexterity and tactile feel. Surgical gloves are available in size 6 thru 8-1/2 and size 9 in two thicknesses.

Developed with the physician in mind, these radiation-attenuating gloves offer the best protection with the greatest dexterity, sensitivity and comfort compared to similar products. Ideal for Radiology, Cath/Electro Physiology, Orthopedic Surgery, Endoscopy, Urology, and Pain Management. Ideal for all high exposure procedures where sensitivity and protection are essential. Color – brown. May be re-sterilized with ethylene oxide gas or Sterad up to three times before disposal. Sizes 6-8-1/2 and 9. Falls under latex-free guidelines.

ATTENUATION—RADIATION REDUCTIONS %				
60 kVp	80 kVp	100 kVp	120 kVp	THICKNESS
43%	35%	26%	23%	0.088"-0.2200mm
55%	43%	35%	31%	0.012"-0.300mm



### AVAILABLE IN 6 SIZES!

COATED		
CAT. NO.	SIZE	PRICE
AX-6.5	6.5	\$32.50
AX-7	7	\$32.50
AX-7.5	7.5	\$32.50
AX-8	8	\$32.50
AX-8.5	8.5	\$32.50
AX-9	9	\$32.50

**TO PLACE AN ORDER, JUST PRINT AND COMPLETE THE FORM BELOW AND FAX TO 302.248.9310 OR SEND TO: [orders@imagingervices.net](mailto:orders@imagingervices.net)**

Please ship me (circle Cat. No.) -

Name/Practice: \_\_\_\_\_

Address: \_\_\_\_\_

CITY/STATE/ZIP: \_\_\_\_\_

CREDIT CARD NO.: \_\_\_\_\_

EXP DATE: \_\_\_\_\_ CODE: \_\_\_\_\_